## DIAVEST TURBO

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#### *Phosphate bound investment for model casting alloys*

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*DIAVEST TURBO is an investment for speed heating procedure*

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*We recommend the use of a suitable silicone duplicating material. Hydro gel duplicating is possible. **200 g powder and 36 ml total mixing liquid** (expansion liquid + water) are required for one model. See schedule overleaf! Apply debubbler (silicone & wetting agent no. 25030) onto the silicone surface. Avoid puddles. The duplicating mould must be absolutely dry when being poured with investment. If necessary, use air pistol softly. Pour investment under slightest possible vibration. Switch vibrator off immediately after pouring. Remove hardened model from the mould after 45 minutes and dry it in a drying unit at 80°C for 40 minutes. Application of an adhesive agent (adhesive spray no. 14140) is recommended for improved adhesion and adaption of wax patterns on the dried out surface of the model. The additional application of a surface hardener is recommended in case of gel duplication.*

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##### **Production of the lining mould and heating sequences**

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*Complete the casting mould / lining with **400 g powder and 72 ml total mixing liquid** (expansion liquid + water). See schedule! Avoid strong and long vibration! Switch vibrator off immediately after pouring. 20 min. after start mixing, store muffle in the preheated furnace at end temperature (950 – 1050°C). After at least 60 min. at end temperature, the muffle can be casted. Let it slowly self-cool down to room temperature. **Do not quench the casting!***

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**Model:** 200 g powder : 27 ml liquid + 9 ml purified H<sub>2</sub>O.  
( 36 ml total mixing liquid )

**Lining:** 400 g powder : 54 ml liquid + 18 ml purified H<sub>2</sub>O.  
( 72 ml total mixing liquid )

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*Mix H<sub>2</sub>O and expansion  
liquid first, then add powder.*

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*Store investment in well-sealed containers. Keep away from moisture.*

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DFS-DIAMON GmbH \* Ländenstraße 1 \* D-93339 Riedenburg \* +49 9442 91890