# BIODURSOFT

CE 0297  
Made in Germany

Ländenstraße 1 | D-93339 Riedenburg  
Phone +49 (0) 94 42 | 9189-0  
info@dfs-diamon.de | www.dfs-diamon.de

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### Usage:

Casting and bonding alloy type 4 for crown- and bridgework for high fusing ceramics, to be used with all qualified, common casting techniques.

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#### 1. Waxup

Wall thickness for single crowns: minimum 0.3 mm  
Wall thickness for bridge crowns: minimum 0.4 mm  
Build an anatomically minimized bridgeframe sufficiently solid and stable relative to its span length. Avoid any sharp edges where porcelain is to be bonded. Incorporate large interdental abutting surface areas, where soldering is planned in spanned bridgework. Use connecting sprues with a diameter of 3.0 mm minimum. Wax T-bars or main bars should have a diameter of 4-5 mm minimum.

#### 2. Wax framework and Investing

Apply „Silicone & Wax Wetting Agent“ (order-no. #25030), or „Waxurit“ (order-no. #17111), in order to ensure a torsion-free removal of wax bridges. Use vacuum stirring unit for mixing phosphate bonded investment (Vesto-Fix). Best results will be achieved when using the „Thermofix 2000“ casting system.

#### 3. Melting and Casting

Only use new material!  
Follow the directions for use of the investment manufacturer.  
Heat the muffle up to 950 °C  
Liquidus point 1395 °C  
Casting temperature 1500 °C  
Only use ceramic casting crucibles. For induction or open flame melting.  
Start casting when the skin of the melting metal starts to tear open.

#### 4. Finishing of the Framework

Do not quench the casting! Cool off slowly to room temperature. Only use carbides like „Diadur“ carbide finishers with medium, fine or superfine crosscut to achieve a clean, contamination-free surface. Maintain speed between 15.000 to 30.000 rpm. Keep driving the carbide finisher in one way motion. **Please note, that to be bonded metal surfaces must not be prepared with sintered diamonds, electroplated diamonds or stones as they will contaminate the alloy with undesirable residue, which cannot be completely removed.** Sandblasting with aluminium oxide (approx. 130 µ) has to be carried out with an impact angle to the metal's surface of 45° in order to guarantee an evenly blasting and to avoid unremovable contamination. Boiling for 10 min. in aqua dest or steamblasting is necessary to successfully complete the finishing. Beware of any contamination of the steam pistol's edge with wax or oil. **Avoid inhaling metal residue!**

#### 5. Degassing the Alloy

Degassing may not be needed when following the preceding steps properly. It can be carried out to control the cleanliness of the metal surface by firing for at least 5 minutes under vacuum at 960 °C. Remove the metal's oxide layer with carbides and aluminium oxide sandblasting completely and finish off again by steaming.

#### 6. Wash Opaque Firing

Fire at 20 °C higher temperature than recommended by the ceramics manufacturer.

#### 7. Opaque Firing

Follow the directions for use of the porcelain manufacturer.

#### 8. Dentine Firing

Follow the directions for use of the porcelain manufacturer. Longtime cooling!

#### 9. Soldering

Keep soldering base as small as possible! Preheat soldering base 10 min. at 500 °C. Follow the directions for use of the solder manufacturer. Use flux and cool to room temperature. Remove conglomerated flux with carbides and aluminium oxide sandblasting.

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