DIADUR

Ländenstraße 1 | D-93339 Riedenburg
Phone +49 (0) 94 42 | 9189-0
info@dfs-diamon.de | www.dfs-diamon.de

CE 0297
Made in Germany

02013

Lega dentale in CoCr CoCr Dental alloy

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Usage:

Nickel free alloy type 5 for partial frame work.

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DIADUR is a high-grade cobalt-chrome alloy for partial frame work.
DIADUR is Ni-free

1. Sprues

Every casting - depending on the dimension of the framework - should be fitted with 2 to 3 sprues: attach wax sprues with 3-4 mm diameter at the thickest points without tapering in order to facilitate the influx and controlled solidification. The length of the sprues must be determined individually by the position of the casting cone.

2. Investment

Release the surface tension at the wax surface with „Silicone and wax wetting agent“ (order-no. #25030). Using the investment „Diavest“ (order-no. #28064-3) will guarantee best results, smooth surfaces and castings with excellent fittings.

3. Preheating (speed)

Follow the instructions of the investment manufacturer.

4. Casting

Do not add more than 20% used alloy to the material to be melted. Cast with any qualified casting unit. Use ceramic crucibles. Start casting before the skin of the melted alloy tears open. Do not use flux.

5. Adjustment and Polishing

We recommend rotaries such as carbides, sinter diamonds, silicone polishers and polishing wheels with short bristles and polishing paste. **Avoid inhaling metal residue!**

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