

FG Diamanten – Royal line in 2-Lagen Diamantierung

FG diamonds – Royal line with two-layer diamond technology

ISO Nr.: 806 315 F 504 superfein / *superfine*

806 315 F 514 fein / *fine*


806 315 F 524 mittel / *medium*

806 315 F 534 grob / *coarse*

806 315 F 544 extragroß / *extracoarse*

Drehzahl / *speed*

30' - 300000 Upm / *rpm*



801L	F :	697	697	697	697	697	697	697
	d :	012	014	016	018	021	023	027
	KL/hl :	1,2	1,4	1,6	1,8	2,1	2,3	2,7
	P :	1	1	1	1	2	2	3

955	F :	699
	d :	008
	KL/hl :	3
	P :	1



FG Diamanten – Royal line in 2-Lagen Diamantierung

mit Tiefenmarkierung zur Inlaypräparation

FG diamonds – Royal line with two-layer diamond technology with depth marker for inlay preparation

ISO Nr.: 806 314 F 504 superfein / *superfine*

806 314 F 514 fein / *fine*


806 314 F 524 mittel / *medium*

806 314 F 534 grob / *coarse*

806 314 F 544 extragroß / *extracoarse*

Drehzahl / *speed*

30' - 300000 Upm / *rpm*



845KR	F :	544
	d :	025
	KL/hl :	4
	P :	2

846KR	F :	545
	d :	018
	KL/hl :	6
	P :	1

959KR	F :	584
	d :	018
	KL/hl :	5,5
	P :	1

532858L

532868L

534128L

Alle Royal Line FG Diamanten sind zu 5 Stück verpackt.
 F = ISO-Figur-Nr. · KL = Kopflänge in mm
 d = Durchmesser in 1/10 mm · P = Preiskategorie

All royal line FG diamonds are packed as 5 pieces.
 F = ISO-Shape-No. · hl = head length in mm
 d = diameter in 1/10 mm · P = price category