




FG Diamanten – Royal line minimal invasiv



FG diamonds – Royal line minimally invasive

ISO Nr.: 806 314 F 514 fein / fine
806 314 F 524 mittel / medium

Drehzahl / speed
30' - 300000 Upm / rpm



		F:	001	001	001	001
		d:	009	012	014	016
		KL/hl:	-	-	-	-
		P:	5	5	2	2
		801				

		F:	002	002	002
		d:	008	010	012
		KL/hl:	2,5	2,5	3
		P:	5	5	5
		802			

537418

537268



537278



537288



537006
537008

537016
537018

537028

		F:	138
		d:	007
		KL/hl:	3
		P:	5
		838	

		F:	156
		d:	008
		KL/hl:	2
		P:	5
		835K	

		F:	237	237	237
		d:	008	010	012
		KL/hl:	2,5	2,5	2,5
		P:	5	5	5
		830			



537318



537106
537108



537116
537118



537126
537128

537398

		F:	238
		d:	009
		KL/hl:	3
		P:	5
		830	

		F:	271
		d:	014
		KL/hl:	3
		P:	2
		890	

		F:	536
		d:	007
		KL/hl:	3
		P:	5
		889	

		F:	540	540
		d:	007	009
		KL/hl:	3,5	3,5
		P:	5	5
		889		

537328

537388

537308

537206
537208

537216
537218

Alle Royal Line FG Diamanten sind zu 5 Stück verpackt.
F = ISO-Figur-Nr. · KL = Kopflänge in mm
d = Durchmesser in 1/10 mm · P = Preiskategorie

All royal line FG diamonds are packed as 5 pieces.
F = ISO-Shape-No. · hl = head length in mm
d = diameter in 1/10 mm · P = price category