

FG Diamanten kurz – Royal line in 2-Lagen Diamantierung


FG diamonds short shank – Royal line with two-layer diamond technology











ISO Nr.: 806 313 F 514 fein / fine








Drehzahl / speed

806 313 F 524 mittel / medium

30' - 300000 Upm / rpm



	  								
		F:	001	001	001	001	001	001	001
		d:	009	010	012	014	016	018	
		KL/hl:	-	-	-	-	-	-	-
P:	1	1	1	1	1	1	1		

	  				
		F:	002	002	002
		d:	009	012	016
		KL/hl:	2,5	3	3
P:	1	1	1		

53528ku

530078ku

530018ku

530028ku






530038ku






530048ku






535648ku






530088ku

530338ku

	 			
		F:	010	010
		d:	010	012
		KL/hl:	1	1,5
P:	1	1		

	 			
		F:	019	019
		d:	010	012
		KL/hl:	1	1,5
P:	1	1		

	 			
		F:	031	031
		d:	031	
		KL/hl:	3,1	
P:	3			

	 			
		F:	108	
		d:	008	
		KL/hl:	3	
P:	1			

530168ku

530108ku

530208ku









530218ku






534366ku






530608ku

New

New

	 						
		F:	109	109	109	109	109
		d:	009	010	012	014	016
		KL/hl:	3	4	4	4	4
P:	1	1	1	1	1		

	 			
		F:	110	110
		d:	012	014
		KL/hl:	6	6
P:	1	1		

	 			
		F:	140	
		d:	013	
		KL/hl:	6,5	
P:	1			

535518ku

530618ku

530615ku

530628ku

530638ku

535678ku

530645ku

530655ku

532463ku

Alle Royal Line FG Diamanten sind zu 5 Stück verpackt.
 F = ISO-Figur-Nr. · KL = Kopflänge in mm
 d = Durchmesser in 1/10 mm · P = Preiskategorie

All royal line FG diamonds are packed as 5 pieces.
 F = ISO-Shape-No. · hl = head length in mm
 d = diameter in 1/10 mm · P = price category